


# **Bacteria Classification**


**Streptococcus Spp.**

**Lab 2**

## Scientific Classification of Streptococcus Spp.

- **Domain: Bacteria**
- **Phylum: Firmicutes**
- **Class: Bacilli**
- **Order: Lactobacillales**
- **Family: Streptococcaceae**
- **Genus: Streptococcus . . . . . include at least 77 species.**





# **Blood Agar**


Several species of Gram-positive cocci produce exotoxins called Hemolysins, which are able to destroy red blood cells (RBCs) and hemoglobin. Blood Agar, sometimes called Sheep Blood Agar because it includes 5% sheep blood in a Tryptic Soy Agar base, allows differentiation of bacteria based on their ability to hemolysis RBCs.

The three major types of hemolysis are  $\beta$ -hemolysis,  $\alpha$ -hemolysis, and  $\gamma$ -hemolysis.

❑  $\alpha$ -hemolysis is the partial destruction of RBCs and produces a greenish discoloration of the agar around the colonies

❑  $\beta$ -hemolysis, the complete destruction of RBCs and hemoglobin, results in a clearing of the medium around the colonies.

❑  $\gamma$ -hemolysis is non hemolysis and appears as simple growth with no change to the medium.



**Hemolysins produced by streptococci are called streptolysins. They come in two forms- type O and type S. Streptolysin**


- **Streptolysin O:** is oxygen-labile and expresses maximal activity under anaerobic conditions.
- **Streptolysin S:** is oxygen-stable but expresses itself optimally under an aerobic conditions as well.

**\*\*\*The easiest method of providing an environment favorable for streptolysins on Blood Agar is what is called the **Streak–stab technique.****


**In this procedure the Blood Agar plate is streaked for isolation and then stabbed with a loop. The stabs encourage **streptolysin activity because of the reduced oxygen concentration of the subsurface environment.****

## The major use of Bile Aesculin Agar

is to differentiate between enterococci / Group D *streptococci* (*Enterococcus*: *E. faecalis* and *E. faecium*) and non-Group D streptococci. It may also be used for the presumptive identification of other groups of organisms.




**Enterococci / Group D streptococci hydrolyse aesculin to form aesculetin and dextrose. Aesculetin combines with ferric citrate in the medium to form a dark brown or black complex which is indicative of a positive result. Bile salts will inhibit Gram-positive bacteria other than enterococci / Group D streptococci. The value of bile tolerance together with hydrolysis of aesculin as a means of presumptively identifying enterococci / Group D streptococci is widely recognized.**





**Streptococcus under light microscope**



**Streptococcus under electronic microscope**



Optochin disk test




Plate on the right is sensitive

Bacitracin test