

Bacteria Classification

Escherichia coli

Lab 8

Scientific Classification of *Escherichia spp.*

Domain: Bacteria

Phylum: Proteobacteria

Class: Gammaproteobacteria

Order: Enterobacteriales


Family: Enterobacteriaceae

Genus: Escherichia

Species: *E. coli*

Enterobacteriaceae

Lactose fermentation on M. A



The IMViC Tests

Indole, Methyl red, Voges-Proskauer, and Citrate; the “i” is for ease of pronunciation.


The amino acid **tryptophan** is found in nearly all proteins. Bacteria that contain the enzyme **tryptophanase** can hydrolyze tryptophan to its metabolic products, namely, indole, pyruvic acid, and ammonia. The bacteria use the pyruvic acid and ammonia to satisfy nutritional needs; indole is not used and accumulates in the medium. The presence of indole can be detected by the addition of **Kovacs' reagent**. Kovacs' reagent reacts with the indole, producing a bright red compound on the surface of the medium. Bacteria producing a red layer following addition of Kovacs' reagent are **indole positive**; the absence of a red color indicates tryptophan was not hydrolyzed, and the bacteria are **indole negative**.




Citrate test



Methyl Red Test



Voges-Proskauer test



Indole Production

Procedure

- a. Label each of the Tryptophan Broth tubes with the name of the bacterium to be inoculated (*E. coli*), (If tryptophan broth medium is unavailable, tryptic soy broth is a good substitute for testing indole production.)
- b. Using aseptic technique, inoculate each tube by a loopful of culture.
- c. Incubate the tubes for about 24 hours at 35°C.
- d. Remove the tubes from the incubator and add **0.5** ml (about **10** drops) of ***Kovacs' reagent*** to each tube, and shake the tube gently. A deep red develops in the presence of indole. Negative reactions remain colorless or light yellow.

2- Methyl Red Test

All enteric bacteria catabolize glucose for their energy needs; however, the end products vary depending on the enzyme pathways present in the bacteria. The pH indicator methyl red detects a pH change to the acid range as a result of acidic end products such as lactic, acetic, and formic acids. This test is of value in distinguishing between *E. coli* (a mixed acid fermenter) and *Enterobacter aerogenes* (a butanediol fermenter).

Mixed acid fermenters such as *E. coli* produce a mixture of fermentation acids and thus acidify the medium.

Butanediol fermenters such as *E. aerogenes* form butanediol, acetoin, and fewer organic acids. The pH of the medium does not fall as low as during mixed acid fermentation. At a pH of 4, the methyl red indicator turns red, a positive methyl red test. At a pH of 6, the indicator turns yellow, a **negative methyl red test**.

Procedure

- a. Label each of the Methyl red - Voges- Proskauer (**MR-VP**) broth- media tubes with the name of the bacterium.

Typical Formula	Gram /litter
Peptone	7.0
Glucose	5.0
Phosphate buffer	5.0
pH 6.9 ± 0.2	



Methyl Red Test

- b. Inoculate each tube with the appropriate bacterium by means of a loop inoculation.**
- c. Incubate all tubes at 35°C for 24 to 48 hours. For slow fermenters, it may take four to five days.**
- d. Transfer $1/3$ of each culture into an empty test tube and set these aside for the Voges-Proskauer test.**
- e. To the $2/3$ of the culture remaining in each tube, add 0.2 ml (about 4 to 5 drops) of methyl-red solution indicator, and read the colour on the surface of the medium immediately, any color change (a red color is MR positive).**

Glucose → Pyruvic acid → Mixed acid fermentation (pH 4.4)
↓
Red color with methyl indicator

3- Voges-Proskauer Test


The **Voges-Proskauer test** identifies bacteria that ferment glucose, leading to **2,3-butanediol** accumulation in the medium. The addition of (40% KOH and a 5% solution of **alpha-naphthol**) will detect the presence of acetoin- a precursor in the synthesis of 2,3-butanediol.

In the presence of the reagents and acetoin, a cherry-red color develops. Development of a red color in the culture medium 15 minutes following the addition of reagent represents a **positive VP test**; absence of a red color is a **negative VP test**.



Procedure

- Use the $\frac{1}{3}$ aliquot from the methyl red test. While wearing disposable gloves, add 0.6 ml of 5% **alpha naphthol** (15 drops) followed by 0.2 mL of 40% **KOH** (5 drops) to each culture.
- Shake vigorously to aerate. Positive reactions occur at once or within 20 minutes and are indicated by the presence of a red color.



Procedure

The **citrate utilization test** determines the ability of bacteria to use **citrate** as a *sole carbon source* for their energy needs. This ability depends on the presence of a **citrate permease** that facilitates transport of citrate into the bacterium. Once inside the bacterium, citrate is converted to pyruvic acid and CO₂. **Simmons citrate agar** slants contain sodium citrate as the carbon source, ammonium (NH₄⁺) as a nitrogen source, and the pH indicator bromothymol blue. This test is done on slants since O₂ is necessary for citrate utilization. When bacteria oxidize citrate, they remove it from the medium and liberate CO₂. CO₂ combines with sodium (supplied by sodium citrate) and water to form sodium carbonate— an alkaline product. This raises the pH, turns the pH indicator to a blue color, and represents a **positive citrate test**; absence of a color change is a **negative citrate test**. Citrate-negative cultures will also show no growth in the medium.

Procedure

- a. Inoculate each bacterium into Simmons citrate agar slants tube by means of a *stab-and-streak* inoculation.
- b. Incubate these cultures for 24 to 48 hours at 35°C.
- c. Examine the slant cultures for the presence or absence of growth and for any change in color from **green to blue**. The development of a deep **blue** color is a positive test.



Citrate test